

Add Hi-Tech Precision to your Surgical Expertise...









• Heated High Pressure Regulator

Technical Specifications of Shalya Endojet - CO₂ Insufflator

Type of Gas	Medical CO ₂				
Pre-setted Modes	Veress Mode, Endoscopy Mode, Pediatric Mode* & Bariatric Mode*				
Set Pressure	Veress Mode: 1-18 mmHg Pediatric Mode: 1-15 mmHg*		Endoscopy Mode: 1-20 mmHg* 1-30 mmHg		
			Bariatric Mode : 1-30 mmHg*		
Gas Flow Rate	Modes	30 LPM 1 to 5 LPM 1 to 10 LPM 1 to 30 LPM		50 LP	M *
	Veress Mode			1 to 5 L	PM
	Pediatric Mode*			1 to 10 LPM	
	Endoscopy Mode			1 to 50 LPM	
	Bariatric Mode*	1 t	o 30 LPM	1 to 50 l	LPM
Flow rate setting	1 LPM/step				
Pressure setting	1 mmHg/step				
Connection Options	Cylinder/Central Gas Supply (Optional)				
Maximum input gas supply pressure	80 Bar / 1160 psi / 60,000 mmHg				
Maximum output gas supply pressure	0.07 Bar / 1.06 psi / 45 mmHg				
Input Voltage	95 to 260 V AC, 50/60 Hz				
Audio Indicator	Speaker to give error tone for safety.				
Touch Control Panel	5" LCD screen				
Design Standard	General Safety Standards: IEC	ral Safety Standards: IEC 60601-1			
	EMI/EMC Compatibility Standard: IEC 60601-1-2				

Features of Shalya Endojet

- **High Flow Insufflation** with maximum insufflation flow rate **upto 30 or 50* LPM**.
- Accurate Flow Control: Flow rate can be set in step of 1 LPM for delicate & critical surgeries.
- Automatic Temperature Control of the CO₂ to avoid Hypothermia & Hyperthermia.
- User Settable Temperature Control for heater system.*
- Clear & Unobstructed View as warmed CO₂ gas avoids fogging over the camera.
- Smoke & Mist Evacuator provides clear & unobstructed view of surgical field during laparoscopic procedures.*
- Surgery Timing Indicator indicates the time from the onset of surgery.
- Over Pressure Detection controls the gas pressure without user's interruption.
- **Gas Indicator** for remnant & consumed gas in the cylinder being displayed on-screen.
- **Cylinder Status** gives continuous visual indication of present state of the cylinder.
- Automatic Pressure Regulation with internal exhaust port which prevents the pneumo-peritonueum from collapsing or overpressure.
- Tube obstruction alarm in case gas tubing gets obstructed or pinched during surgery.*
- **Virtual/External Pressure Alarm**: Unit gives an audio visual alarm & holds the abdominal pressure when sudden external pressure is applied during trocar insertion.
- Surgery History maintains upto last 100 surgery records for Medico legal compliance.
- Control via USB/Ethernet/CAN for centralized control in conjugation with other Shalya OT equipments.

Applications of Shalya Endojet

Widely used for Surgical Laparoscopic procedures like Total Laparoscopic Hysterectomy (TLH), Lap Cholecystectomy, Myomectomy, Lap Appendectomy, Laparoscopic Hernia Repair, Laparoscopic Hepatectomy, Laparoscopic Pancreatectomy, Laparoscopic General Surgery & Laparoscopic Cancer Surgery.









Note - * Features not applicable for Endojet E

XcelLance Medical Technologies Pvt. Ltd.

Available at: Xavier Med Pvt. Ltd.

1-10-51/2, Road 6B, Dwarakapuram, Dilsukhnagar, Hyderabad - 500060

Customer Care: 8008895680, Service: 8008895079