## **Basic Anaesthesia Machine**





## Features:

- Trolley mode of M.S. power coated, stainless steel table, top tray, modified monitoring accessories tray, drawer top stainless tray.
- Anti Hypoxic Rota meter its pneumatic device ensure no nitrous oxide flow, if oxygen flow is lower than 0.5 ltr/min approx.
- Minimum 25% (+/-5%) oxygen in surgery gas mixture.
- Ratio Control 1:3 mixed ratio between 0<sub>2</sub>& NQ gas flow.
- · Oxygen driven NO supply.
- Oxy failure warning device (OFWD).
- · Higher efficiency forged regulators & yokes with SS fittings.
- · High precision pressure relief cum non return valve
- · Vaporizers Gold man Halothane.
- Change over mechanism with oxygen flush (Open & close circuit)
- Gas inlet Port O<sub>2</sub>& NQ each 3 nos. Pin-indexed gas cylinder Connections for O, & NO each 2 nos & diameter index gas specific central pipe line connection each 1 no for O<sub>2</sub>& NQ
- Standard Magills circuit, all rubber antistatic tested. 5" dia antistatic breaking castor wheel.
- Oxygen outlet for driving ventilator (Double).
- 10 ltr W/C cylinder carrier facility, Provision for Tec Vaporizer
- Circle Absorber & Ventilator, (Optional at extra cost)
- Size 1490 mm x 910 mm x 750 mm (HxWxD)

Available at:

## XAVIER MED PVT. LTd.

1-10-51/2, Road 6B, Dwarakapuram, Dilsukhnagar, Hyderabad - 500060 www.buyxayier.com

Customer Care: 8008895680, Service: 8008895079