

# Basic Anaesthesia Machine



Embedded  
Pressure Gauges

Model:  
**COMPACT**

## Features:

- Anti Hypoxic pneumatic device. Ensure no nitrous oxide flow, if oxygen flow is lower than 0.5 Ltr / min approx.
- Minimum 25% (+/-5%) oxygen in surgery gas mixture
- Ration control 1:3 mixed ratio between O<sub>2</sub>, & N<sub>2</sub>O gas flow.
- Rota meter 9"
- Oxygen driven N<sub>2</sub>O supply.
- Audio visual oxy failure warning device (OFWD)
- Higher efficiency forged regulators & yokes with SS fittings.
- High precision pressure relief cum non return valve
- Vaporizer Gold Man Halothane
- Change over mechanism with oxygen flush. (open & close circuit)
- Gas inlet Port O<sub>2</sub> & N<sub>2</sub>O each 3 Nos
- Pin-indexed gas cylinder connections for O<sub>2</sub> & N<sub>2</sub>O each 2 Nos & diameter index gas specific central pipe line connection each one no. for O<sub>2</sub> & N<sub>2</sub>O
- Standard Mogills circuit all rubber antistatic tested.
- Trolley mode of M.S multi colour power coated.
- Stainless steel table top tray
- Modified monitoring accessories tray, drawer top stainless steel tray
- 5" dia antistatic breaking castor wheel
- Oxygen outlet for driving ventilator (double)
- Large dia (63mm) panel mounted cylinder pressure, operating pressure & pipeline pressure display gauges for O<sub>2</sub> & N<sub>2</sub>O
- 10 Ltr. W/C cylinder carrier facility. Provision for tec vaporizer, circle absorber & ventilator. (optional at extra cost)
- Size 1530 mm x 900 mm x 770 mm (H x W x D)

Available at:  
**XAVIER Med Pvt. Ltd.**  
1-10-51/2, Road 6B, Dwarakapuram, Dilukhnagar, Hyderabad - 500060  
[www.buyxavier.com](http://www.buyxavier.com)  
Customer Care: 8008895680, Service: 8008895079