

HiOxy 1570s

High Flow Oxygen Therapy System

HiOxy 1570s High Flow oxygen therapy is a simple to use system that delivers warmed and humidified respiratory gases to spontaneously breathing patients and is capable of providing High Flow Therapy as well as Bubble CPAP.



Dual flow blender



Nasal cannula (Neonate)



Nasal cannula (Pediatric)



Nasal cannula (Adult)

FOG 2140A Servo Controlled Humidifier

(Heater with sensor cable & adapter included)

- Suitable for: Adults, Children.
- Mode of application: Invasive and noninvasive.
- Display: 2.8" high resolution TFT color screen display.
- Rated voltage: **220VAC/50Hz.**
- **Heating plate power: 150W.**
- **Heating coil power:**
 - A. The inspiration side: 24VDC, 36W;**
 - B. The expiration side: 24VDC, 36W;**
- * setting range for heating temperature (adjustment accuracy of 0.5 °C)
 - A. Invasive mode: patient mode end temperature: 35-40 °C, the**
 - B. Noninvasive mode: patient mode end temperature: 30-37 °C,**
- Temperature display error: ± 0.5 °C
(temperature in range of 22-45 °C, under stable airflow).
- Humidity control range: $\geq 95\%$ RH.
- Machine noise: ≤ 50 Db.



Humidifier / Temperature sensor feedback cable



Blender

It features a 0-70 lpm flow meter on the left side and a 0-15 lpm flow meter on the right side.



SPECIFICATION

Oxygen Concentration Range	21-100%
Number of Ports	3 (2 occupied by flow meters plus a bottom port)
Left Flowmeter Range	15-70lpm (2-70 if bleed flow is on right flow meter is in upright position)
Right Flowmeter range	2-15 LPM
Accuracy	$\pm 3\%$ Full scale
Weight	3.17lbs (1.44kg)
Size	7.98W x 6.64H x 4.56D in (20.3 x 16.9 x 11.6 cm)



Av. Gal. Ataliba Leonel, 1790
CEP 02033-020 - São Paulo - SP - Brasil
PABX: 55 (11) 2972-5700 - FAX: 55 (11) 2979-1575
Email: comercial@fanem.com.br - Site: www.fanem.com.br

Available at:

XAVIER Med Pvt. Ltd.

1-10-51/2, Road 6B, Dwarakapuram,
Dilsukhnagar, Hyderabad - 500060

www.buyxavier.com

Customer Care: 8008895680,
Service: 8008895079