

# BiPAP - BiLEVEL ST



Modes	CPAP / AutoCPAP / AutoBiLEVEL / BiLEVEL S / BiLEVEL T / BiLEVEL ST
Working Pressure Range	4-35 mBar
SD Memory Card	4 GB
Usage Features	3.5 inch Screen Auto FRT (Flexible Respaps® Therapy) Ramp (0-60 min) Mask Test (0-45 Sec) HB (Fast Start) Tubing Selection Mask Selection Wake Up Alarm Automatic Accessory Detection Screen Timeout
Additional Indicators	Pressure Frequency (4-40 BPM) I/E Ratio Leakage Ratio PF (Peak Flow) Volume (200-1500ml) MV (Minute Ventilation) Graphical Screen (Real Time Data Tracking)
Humidifier	Embedded Humidifier 10 Level Heating
Alarms	Apnea Alarm & Leakage Alarm
Menus	Settings Change Screen Reports Info Alarm History

## BIPAP / CPAP Circuit

Optional Hospital Special Circuit for Prevention of Patient Exhaled contaminated Air into Ward.



Including Mask with Oxygen Enrichment Port, Tubing with Bacterial / Viral Filter, Catheter Mount with Y Connector & Bacterial Viral Filter.

Available at:

**XAVIER Med Pvt. Ltd.**  
1-10-51/2, Road 6B, Dwarakapuram,  
Dilsukhnagar, Hyderabad - 500060  
[www.buyxavier.com](http://www.buyxavier.com)  
Customer Care: 8008895680,  
Service: 8008895079