

CPAP

with Flexible Respaps® Therapy



Modes	CPAP
Working Pressure Range	4-20 mBar
SD Memory Card	4 GB
Usage Features	3.5 inch Screen Auto FRT (Flexible Respaps® Therapy) Ramp (0-60 min) Mask Test (0-45 Sec) HB (Fast Start) Tubing Selection Mask Selection Wake Up Alarm Automatic Accessory Detection Screen Timeout
Additional Indicators	Pressure Graphical Screen (Real Time Data Tracking)
Humidifier	Embedded Humidifier 10 Level Heating
Alarms	Leakage Alarm
Menus	Settings Change Screen Reports Info Alarm History

CPAP / BIPAP Circuit

Optional Hospital Special Circuit for Prevention of Patient Exhaled contaminated Air into Ward.



Including Mask with Oxygen Enrichment Port, Tubing with Bacterial / Viral Filter, Cather Mount with Y Connector & Bacterial Viral Filter.

Available at:

XAVIER Med Pvt. Ltd.
1-10-51/2, Road 6B, Dwarakapuram,
Dilsukhnagar, Hyderabad - 500060
www.buyxavier.com
Customer Care: 8008895680,
Service: 8008895079